

# HIPOKRAT 2000

KNEE REPLACEMENT SYSTEM

## CRUCIATE RETAINING AND POSTERIOR STABILIZED OPTIONS

9 sizes (right / left), CoCrMo alloy

Deep patella groove provides stable patellofemoral tracking

Rough surface provides strong cementation

5 sizes, 6 thickness  
(6, 8, 10, 12, 15, 18 mm)

9 sizes metal back offers optimal fit

Higher anterior surface will prevent the paradoxical anterior sliding

Cemented, uncemented and MRI safe conditional options