



CEMENTED MODULAR STRAIGHT STEM
ÇİMENTOLU MODÜLER DÜZ STEMMaterial: CoCrMo
Cone: 12/14

Code Kod	Size Boy
010.3600	7.5 mm
010.3700	10 mm
010.4000	11.25 mm
010.3800	12.5 mm
010.4100	13.75 mm
010.3900	15 mm
010.4200	16.25 mm
010.4300	17.5 mm
010.4400	20 mm



Hipokrat Tıbbi Malzemeler İmalat ve Pazarlama A.Ş.

📍 7407/6 Sokak No:10 35060 Pınarbaşı Bornova - İZMİR / TÜRKİYE

☎ + (90) 232 479 11 09

✉ info@hipokrat.com.tr

☎ + (90) 232 479 44 29

🌐 www.hipokrat.com.tr