



Çimentosuz Asetabular Vidalı Kap
Uncemented Press-fit Acetabular Cup

Material: Titanium Alloy



TIPS+HA Coated
TIPS+HA Kaplamalı

TIPS+HA TIPS+HA	Diameter Çap	Linear Size Linear Çap
077.1600	ø42mm	AA
077.1700	ø44mm	AA
077.1800	ø46mm	A
077.1900	ø48mm	A
077.2000	ø50mm	B
077.2100	ø52mm	B
077.2200	ø54mm	C
077.2300	ø56mm	C
077.2400	ø58mm	D
077.2500	ø60mm	D
077.2600	ø62mm	E
077.2700	ø64mm	E
077.2800	ø66mm	F
077.2805	ø68mm	F

Material: Titanium Alloy



TIPS Coated
TIPS Kaplamalı

TIPS TIPS	Diameter Çap	Linear Size Linear Çap
077.0100	ø42mm	AA
077.0200	ø44mm	AA
077.0300	ø46mm	A
077.0400	ø48mm	A
077.0500	ø50mm	B
077.0600	ø52mm	B
077.0700	ø54mm	C
077.0800	ø56mm	C
077.0900	ø58mm	D
077.1000	ø60mm	D
077.1100	ø62mm	E
077.1200	ø64mm	E
077.1300	ø66mm	F
077.1305	ø68mm	F

Material: Titanium Alloy



Sand Blasted
Kumlamalı

Sand Blasted Kumlama	Diameter Çap	Linear Size Linear Çap
077.0405	ø42mm	AA
077.0410	ø44mm	AA
077.0415	ø46mm	A
077.0420	ø48mm	A
077.0425	ø50mm	B
077.0430	ø52mm	B
077.0435	ø54mm	C
077.0440	ø56mm	C
077.0445	ø58mm	D
077.0450	ø60mm	D
077.0455	ø62mm	E
077.0460	ø64mm	E
077.0465	ø66mm	F
077.0470	ø68mm	F

Material: Titanium Alloy



Multi Holes TIPS+HA Coated
Çoklu Delik TIPS+HA Kaplamalı

TIPS+HA TIPS+HA	Diameter Çap	Linear Size Linear Çap
077.1600	ø42mm	AA
077.1700	ø44mm	AA
077.1800	ø46mm	A
077.1900	ø48mm	A
077.2000	ø50mm	B
077.2100	ø52mm	B
077.2200	ø54mm	C
077.2300	ø56mm	C
077.2400	ø58mm	D
077.2500	ø60mm	D
077.2600	ø62mm	E
077.2700	ø64mm	E
077.2800	ø66mm	F
077.2805	ø68mm	F

Material: Titanium Alloy



Multi Holes TIPS Coated
Çoklu Delik TIPS Kaplamalı

TIPS TIPS	Diameter Çap	Linear Size Linear Çap
077.0100	ø42mm	AA
077.0200	ø44mm	AA
077.0300	ø46mm	A
077.0400	ø48mm	A
077.0500	ø50mm	B
077.0600	ø52mm	B
077.0700	ø54mm	C
077.0800	ø56mm	C
077.0900	ø58mm	D
077.1000	ø60mm	D
077.1100	ø62mm	E
077.1200	ø64mm	E
077.1300	ø66mm	F
077.1305	ø68mm	F

Material: HIGH CORSS LINKED Polyethylene



Liner for Acetbular Cup
Asetabular Kap Yatağı

0° 0°	10° 10°	20° 20°	Linear Size Linear Çap	
101.0837	101.0700	102.0100	AA	ø42-44 / 26 mm
101.0839	101.0102	102.0102	A	ø46-48 / 26 mm
101.0831	101.0100	102.0200	A	ø46-48 / 28 mm
101.0832	101.0200	102.0300	B	ø50-52 / 28 mm
101.0833	101.0300	102.0400	C	ø54-56 / 28 mm
101.0834	101.0400	102.0500	D	ø58-60 / 28 mm
101.0835	101.0500	102.0600	E	ø62-64 / 28 mm
101.0836	101.0600	102.0700	F	ø66-68 / 28 mm
101.0841	101.0801	101.0811	B	ø50-52 / 32 mm
101.0842	101.0802	101.0812	C	ø54-56 / 32 mm
101.0843	101.0803	101.0813	D	ø58-60 / 32 mm
101.0844	101.0804	101.0814	E	ø62-64 / 32 mm
101.0845	101.0805	101.0815	F	ø66-68 / 32 mm
101.0851	101.0821	101.0825	D	ø58-60 / 36 mm
101.0852	101.0822	101.0826	E	ø62-64 / 36 mm
101.0853	101.0823	101.0827	F	ø66-68 / 36 mm

Material: UHMW-PE Polyethylene

0° 0°	10° 10°	20° 20°	Linear Size Linear Çap	
094.1721	094.0700	094.0900	AA	ø42-44 / 26 mm
094.1745	094.0102	094.1102	A	ø46-48 / 26 mm
094.1722	094.0100	094.1100	A	ø46-48 / 28 mm
094.1723	094.0200	094.1200	B	ø50-52 / 28 mm
094.1724	094.0300	094.1300	C	ø54-56 / 28 mm
094.1725	094.0400	094.1400	D	ø58-60 / 28 mm
094.1726	094.0500	094.1500	E	ø62-64 / 28 mm
094.1727	094.0600	094.1600	F	ø66-68 / 28 mm
094.1731	094.1701	094.1711	B	ø50-52 / 32 mm
094.1732	094.1702	094.1712	C	ø54-56 / 32 mm
094.1733	094.1703	094.1713	D	ø58-60 / 32 mm
094.1734	094.1704	094.1714	E	ø62-64 / 32 mm
094.1735	094.1705	094.1715	F	ø66-68 / 32 mm
094.1741	101.0861	101.0865	D	ø58-60 / 36 mm
094.1742	101.0862	101.0866	E	ø62-64 / 36 mm
094.1743	101.0863	101.0867	F	ø66-68 / 36 mm

Material: Titanium Alloy



6,5 mm Titanium Screw
6.5 mm Titanyum Vida

Code Kod	Size Boy
091.1300	ø6.5 x 15mm
091.1400	ø6.5 x 20mm
091.1500	ø6.5 x 25mm
091.1600	ø6.5 x 30mm
091.1700	ø6.5 x 35mm
091.1800	ø6.5 x 40mm
091.1900	ø6.5 x 45mm
091.2000	ø6.5 x 50mm
091.2100	ø6.5 x 55mm
091.2200	ø6.5 x 60mm
091.2300	ø6.5 x 65mm



Hipokrat Tıbbi Malzemeler İmalat ve Pazarlama A.Ş.

📍 7407/6 Sokak No:10 35060 Pınarbaşı Bornova - İZMİR / TÜRKİYE

☎ + (90) 232 479 11 09

✉ info@hipokrat.com.tr

☎ + (90) 232 479 44 29

🌐 www.hipokrat.com.tr