



HBRS Tibia Üst Uç Rezeksiyon Protezi
HBRS Proximal Tibia Resection Prosthesis



**Distal Femur Minimal
Resection Component**
Femur Alt Uç Minimal
Rezeksiyon Parçası

Material: Co. Cr. Mo. Alloy

Left Sol	Right Sağ	Size Boy
647.0409	647.0410	Pediatric
647.0401	647.0402	1
647.0403	647.0404	2
647.0405	647.0406	3
647.0407	647.0408	4



Yoke
Mafsal

Material: Co. Cr. Mo. Alloy

Short Kısa	Long Uzun	Size Boy
647.0908	647.0909	Pediatric
647.0901	647.0902	1
647.0903	647.0904	2



Boss Core Bar
Göbek Mil

Material: Co. Cr. Mo. Alloy

Code Kod	Size Boy
647.0306	28 mm.
647.0303	32,7 mm.
647.0304	36,7 mm.



Bolt for Yoke
Mafsal Vidası

Material: Titanium Alloy

Code Kod
647.0906



Bolt for Small Conical Stem
Stem Vidası

Material: Titanium Alloy

Code Kod
647.0501



Tube for Yoke
Mafsal Tüpü

Material: UHMWPE

Code Kod	Size Boy
647.0911	Pediatric
647.0905	1



Insert
Tibia Yatağı

Material: UHMWPE

6 mm.	10 mm.	Size Boy
647.0605	647.0606	Pediatric
647.0601	647.0602	1
647.0603	647.0604	2
647.0609	647.0610	3
647.0611	647.0612	4



Bush
Burç

Material: UHMWPE

Code Kod	Size Boy
647.0204	Pediatric
647.0202	1



Proximal Tibia Full Resection Component
Proksimal Tibia Ful Rezeksiyon Komponent

Material: Co. Cr. Mo. Alloy

Code Kod	Size Boy
647.0802	Pediatric
647.0801	1



Extension Piece
Uzatma Borusu

Material: Co. Cr. Mo. Alloy

Code Kod	Length Uzunluk
647.1151	40 mm.
647.1152	50 mm.
647.1153	60 mm.
647.1154	70 mm.
647.1155	80 mm.
647.1156	90 mm.
647.1157	100 mm.
647.1158	110 mm.
647.1159	120 mm.

Material: Co. Cr. Mo. Alloy



**Cemented Small
Conical Stem**
Çimentolu Küçük
Konik Stem

Straigh Düz	Ø
647.1101	Ø10 mm.
647.1102	Ø11 mm.
647.1103	Ø12 mm.
647.1104	Ø13 mm.
647.1105	Ø14 mm.
647.1106	Ø15 mm.
647.1107	Ø16 mm.

Material: Titanium Alloy
TIPS Coated



Straigh Düz	Curved Eğri	Ø
647.1110	647.1130	Ø10 mm.
647.1111	647.1131	Ø11 mm.
647.1112	647.1132	Ø12 mm.
647.1113	647.1133	Ø13 mm.
647.1114	647.1134	Ø14 mm.
647.1115	647.1135	Ø15 mm.
647.1116	647.1136	Ø16 mm.

Material: Titanium Alloy
TIPS + HA Coated



Straigh Düz	Curved Eğri	Ø
647.1120	647.1140	Ø10 mm.
647.1121	647.1141	Ø11 mm.
647.1122	647.1142	Ø12 mm.
647.1123	647.1143	Ø13 mm.
647.1124	647.1144	Ø14 mm.
647.1125	647.1145	Ø15 mm.
647.1126	647.1146	Ø16 mm.

**Uncemented Small
Conical Stem**
Çimentosuz Küçük
Konik Stem

Material: Co. Cr. Mo. Alloy



Cemented Large Conical Stem
Çimentolu Büyük Konik Stem

Straight Düz	Curved Eğri	Ø
647.1010	647.1001	Ø10 mm.
647.1011	647.1002	Ø11 mm.
647.1012	647.1003	Ø12 mm.
647.1013	647.1004	Ø13 mm.
647.1014	647.1005	Ø14 mm.
647.1015	647.1006	Ø15 mm.
647.1016	647.1007	Ø16 mm.

Material: Titanium Alloy
TIPS Coated



Straight Düz	Curved Eğri	Ø
647.1030	647.1020	Ø10 mm.
647.1031	647.1021	Ø11 mm.
647.1032	647.1022	Ø12 mm.
647.1033	647.1023	Ø13 mm.
647.1034	647.1024	Ø14 mm.
647.1035	647.1025	Ø15 mm.
647.1036	647.1026	Ø16 mm.

Material: Titanium Alloy
TIPS+HA Coated



**Uncemented Hexagonal
Large Conical Stem**
**Çimentosuz Altıgen
Büyük Konik Stem**

Straight Düz	Curved Eğri	Ø
647.1050	647.1040	Ø10 mm.
647.1051	647.1041	Ø11 mm.
647.1052	647.1042	Ø12 mm.
647.1053	647.1043	Ø13 mm.
647.1054	647.1044	Ø14 mm.
647.1055	647.1045	Ø15 mm.
647.1056	647.1046	Ø16 mm.



Hipokrat Tıbbi Malzemeler İmalat ve Pazarlama A.Ş.

📍 7407/6 Sokak No:10 35060 Pınarbaşı Bornova - İZMİR / TÜRKİYE

☎ + (90) 232 479 11 09

✉ info@hipokrat.com.tr

☎ + (90) 232 479 44 29

🌐 www.hipokrat.com.tr